

## 3.6 Concentration of solutions

**Concentration** is the amount of solute dissolved in a unit volume of solution. The amount of solute is usually measured in moles and the unit volume is usually one litre ( $1 \text{ dm}^3$ ).

The mass of one **mole** of sodium chloride,  $\text{NaCl}$ , is  $58.5 \text{ g}$ . This is the **molar mass** of sodium chloride. To prepare a sodium chloride solution with concentration of exactly  $1 \text{ mol dm}^{-3}$ , one mole (that is  $58.5 \text{ g}$ ) of sodium chloride is placed in a volumetric flask and distilled water added until the volume of the solution is exactly  $1 \text{ dm}^3$ .

### Aim


To prepare  $1 \text{ dm}^3$  of a sodium chloride solution of concentration  $1 \text{ mol dm}^{-3}$ .

### Apparatus and materials

Volumetric flask, $1 \text{ dm}^3$	Funnel
Stopper	Glass rod
Beaker, $100 \text{ cm}^3$	Dropping pipette
Balance	Sodium chloride, $58.5 \text{ g}$
Wash bottle	Distilled water

### Procedure

- 1 Weigh out accurately into a beaker  $58.5 \text{ g}$  of sodium chloride.
- 2 Add the sodium chloride to about  $50 \text{ cm}^3$  of distilled water. Dissolve the solute by stirring with the glass rod.
- 3 Pour the sodium chloride solution into the volumetric flask using a funnel. Add distilled water to the undissolved solid in the beaker and pour this into the flask too.
- 4 Rinse out the beaker with more distilled water and pour the washings through the funnel into the flask.
- 5 Add more distilled water to the volumetric flask. Use a dropping pipette to add the final few  $\text{cm}^3$  to make up the solution to the  $1 \text{ dm}^3$  mark on the flask.
- 6 Stopper the flask and shake the solution.



### Extra Work

- Make magnesium sulphate and potassium chloride solutions with concentrations of  $1 \text{ mol dm}^{-3}$ .
- Make sodium chloride solutions of concentration  $1 \text{ mol dm}^{-3}$  using  $100 \text{ cm}^3$  and  $250 \text{ cm}^3$  volumetric flasks and the appropriate amounts of solute.

### Questions

- 1 a) Why is the distilled water not heated to make the sodium chloride dissolve faster?
  - b) Why is a volumetric flask washed out with distilled water before use?
  - c) Why is distilled water used in this experiment rather than tap water?
  - d) Why is it necessary to shake the volumetric flask at the end of the preparation?
  - e) Why is a dropping pipette, and not the wash bottle, used to make the final few  $\text{cm}^3$  of solution up to the mark on the volumetric flask?
- 2 Copy and complete the following table by calculating the mass of the given compound needed to make a solution of the volume and concentration shown.

compound	volume of solution/ $\text{dm}^3$	concentration of solution/ $\text{mol dm}^{-3}$	mass/g
$\text{Na}_2\text{CO}_3$	1.00	0.50	
$\text{AgNO}_3$	2.00	0.10	
$\text{NaOH}$	0.50	2.00	
$\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$	2.00	0.20	
$\text{KCl}$	5.00	0.25	

- 3 What is the concentration of the following solutions in  $\text{mol dm}^{-3}$ ?

- a)  $1.60 \text{ g}$  of  $\text{Na}_2\text{CO}_3$  in  $200 \text{ cm}^3$  of solution.
- b)  $2.00 \text{ g}$  of  $\text{NaOH}$  in  $1 \text{ dm}^3$  of solution.
- c)  $36.50 \text{ g}$  of  $\text{HCl}$  in  $500 \text{ cm}^3$  of solution.
- d)  $3.32 \text{ g}$  of  $\text{KI}$  in  $100 \text{ cm}^3$  of solution.
- e)  $7.90 \text{ g}$  of  $\text{HNO}_3$  in  $250 \text{ cm}^3$  of solution.